STERYBOX[®] *mini*

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- Medical Environments
- Dental Offices
- Veterinary
- Diagnostic Centers
- Microbiological Labs
- Ambulance
- Nurseries
- Home Care Treatment
- Waiting Rooms
- Quarantine Areas
- Aesthetic & Fitness Centers
- Food Industry
- Pharmaceutical Industry
- Safe Packaging

Effective
Against
SARS virus

progettata e brevettata da
engineered and patented by

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STERYBOX is an innovative device for reducing the bacterial charge of the air in environments at risk of aerial contagion, such as sterilization rooms, dentists' surgeries, laboratories for diagnostics and analysis, operating theatres, first aids, etc. This device uses the well-known properties of the UV-C germicidal radiation at 254 nanometres (UVGI) and, specifically, the capacity of denaturing the genetic material of the cells it strikes, by inactivating every type of microorganism present in the air (bacteria, germs, moulds, viruses). The germicidal radiation, which is isolated inside the device, prevents the spreading of the dangerous UV-C rays in the environment and it eliminates any side effect detrimental to health; the air is aspirated by the device through noiseless axial fans and it is then submitted to an intense UV-C radiations field through a patented coil system; the air is purified and ready to be released into the environment, thanks to a continuous change of air.

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MAIN FEATURES


- Reduces within 40 minutes the aerial contamination into a room up to 99,99%*
- Scientific tests performed by independent laboratories show its efficacy against all viruses, bacteria, spores*
- The worldwide patented system (UV-CHRIS) increases efficacy up to 400% vs all other devices on the market*
- It can work 24 hours a day with in presence of people: no fraction of the germicidal radiation can in fact strike directly eyes and skin
- Specially designed UVGI lamps do not create ozone for a total operative safety*
- On/Off cycles completely automatic, easy installation and safe maintenance, very noiseless even in small rooms*
- Reduces use of surface disinfectants in the office, does not nebulize chemicals in the air, no filters to change
- Installed in several hospitals and surgical theaters in respect of the international regulations (CE marked)

*Tests available upon request

99,99%
Germicidal Reduction
in a Room up to 50 m³
in only 40 minutes


➤ Device for airborne disinfection in medical environments, hospitals, surgeries, dental offices, sterilisation and diagnosis rooms.

Patent system UV-CHRIS (UltraViolet High Reflecting Coil Irradiation System)

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The UV-CHRIS - UltraViolet High Reflecting Coil Irradiation System is a new, patented coil system with a concentration of ultraviolet radiations up to four times higher compared with the standard linear flow devices available on the market. Thanks to it a 50 watt irradiation power is applied on a generation of a linear length of more than 1 metres, with an exposure of the air in direct contact of the lamps that is able to eliminate up to 99.99%* (tests show a virtual 100%) of bacteria and viruses present in the air. This system combines optimal performances with an incredibly small size: width 30 cm, height 70 cm, thickness 15 cm, allowing its installation also in dental offices and sterilization rooms.

Escherichia Coli 8.000
UFC/mc
su M.H. Agar senza
Sterybox
(ricerca lab. Sigla srl)

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Escherichia Coli 8.000
UFC/mc
on M.H. Agar with
Sterybox:
reduction of bacteria is
100%

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Major benefits:

- 1) It can be kept working 24 hours a day regardless of the presence of persons: no fraction of the germicidal radiation can in fact strike directly eyes and skin, the radiation being contained inside the device; unlike the systems with UV-C radiations towards the ceiling, this device guarantees a total operative safety for both the staff and the patients;
- 2) One noiseless axial fans Service-Free of 120 mm with real capacity of 100 cubic metres per hour (capacity in free air more than 150 mc/h), guarantee the airflow for several changes of air. To insure a correct change of air every hour, a device every 50 mc of volume should be installed (a 215 ftc room), multiple devices may be installed: up to 50 mc one device, from 50 to 100 mc two devices on two opposite sides, from 100 to 150 mc three devices on three sides, etc.
The device only needs electricity supply and it can be installed on the wall in few minutes. A silencer reduces the noise made by air coming out of the device, which is lower than an ordinary air conditioning system, so that it is not annoying for the medical staff even in case of continuous working conditions;
- 3) A practical control system with a lead hour-counter ensures a correct maintenance program and consequently an optimal functioning over time. Maintenance is minimum: a change of lamps every 8,000 hours, with an average duration of 2 years for a day use and of 1 year for a continuous use day/night. Lamps are easily replaced in few minutes and original spare parts are labelled with a warranty label (the new lamps cannot be distinguished at sight from the used ones);

www.STERYBOX.com