

DiosolGenerator™

MF Multifunction

For room, air and surface disinfection

DIOP GmbH & Co KG
Disinfection is our Passion



The DiosolGenerator™ disinfects complete environments incl. hard to reach gaps, corners and edges!

Successfully used for...

Disinfection of hospitals and clinics, ambulances, livestock farming, retirement homes, dental practices, building restoration, pharmaceutical industry, food production....

A germ reduction up to 7-log steps is possible!

Low application costs and significant time savings for final disinfection! Low consumption of Diosol disinfectants

Disinfection procedure is characterized by high material compatibility, including sensitive electronics!

The device's features such as mobility and ergonomic design, allow a user-friendly disinfection.

Based on standardized test germs, the disinfection method allows an easy and quick validation.

The spectrum of activity (ABCD) ranges from bacteria, fungi, yeasts, to enveloped viruses and non-enveloped viruses up to spores.

Our products have been tested according to European standards.



Dimensions (W x L x H):	215 x 345 x 330 mm
Weight:	15 kg (excluding disinfectant)
Voltage:	230 V
Consumption:	350 W
Compressor:	Dry air compressor, without oil
Room volume:	270 m ³ (Upgrade up to 450 m ³ and 850 m ³)
Disinfectants consumption:	4 ml / m ³
Noise level:	62 dB

Стерамед ООД

София, ул. Христо Ценов №4
т.: 02 437 24 85
e: office@steramed.bg
i: http://steramed.bg
Дистрибутор за България

DIOP GmbH & Co. KG
Dieselstr. 5-6
D – 61191 Rosbach / Germany
Phone.: +49(0)6003 / 814-800
Fax: +49(0)6003- / 814-671
E-Mail: info@diopgmbh.com
Internet: www.diopgmbh.com