

## МЕДИЦИНСКИ РЪКАВИЦИ – НЕСТЕРИЛНИ – НИТРИЛ



### BIOSOFT PLUS

DAR NI

#### Ръкавици – Нитрил, без талк

Нитрилови ръкавици без талк, светло син цвят, с микро-текстура на пръстите.

**ПРЕДИМСТВА:** висока степен на защита и комфорт, безопасно захващане

- Нитрил, без талк, нестерилни
- Качество - AQL 1.0
- 240 mm маншет
- Светло син цвят
- Микро-текстура на пръстите (грапавост)
- Подходящи за лява и за дясна ръка
- Не пропускащ маншет
- Тествани по: EN 455 I-II-III, EN 374 I-II-III, ASTM F 1670-1671
- DM Class I Dir. 93/42/CEE Dir. 2007/47/EEC, DPI IIIa Cat. CE0465 Dir. EEC 89/686
- Размери: S 6-6,5 / M: 7-7,5 / L: 8-8,5 / XL: 9-9,5
- Опаковка: кутия - дозатор 100 ръкавици - кашон 1,000 ръкавици



### BIOSOFT PLUS

DAR NI

Nitrile glove, light blue colour, powder free, finger texturisation.

#### BENEFITS:

High protection, safe grip and highest comfort.

- Nitrile, Powder free, Non sterile
- AQL 1.0
- 240 mm cuff
- Internally chlorinated
- Light blue colour
- Finger texturisation
- Ambidextrous
- Anti-tear cuff
- DM Class I Dir. 93/42/CEE Dir. 2007/47/EEC, DPI IIIa Cat. CE0465 Dir. EEC 89/686
- Sizes XS: 5-5,5 / S: 6-6,5 / M: 7-7,5 / L: 8-8,5 / XL: 9-9,5
- Dispenser 100 gloves - Carton 1000 gloves