

ХИРУРГИЧНИ РЪКАВИЦИ –СТЕРИЛНИ – ЛАТЕКС, С ТАЛК



FEELTOUCH 300PP SGA1SW

Natural latex Surgical glove, internal coating with modified cornstarch (USP) that is used as a lubricant for easy fit.

Feeltouch offers an excellent combination of protection and comfort, since fits the hand perfectly thanks to its elasticity.

- Latex
- AQL 1,0
- Powdered
- Natural colour
- Micro textured surface
- Anatomic shape (L/R) with curved fingers
- Sterilized by gamma irradiation
- Anti-tear cuff
- Medical Device Class IIa 93/42/EEC
- Test: EN 455 I-II-III, ASTM F 1671
- Sizes 5,5-6-6,5-7-7,5-8-8,5-9
- Dispenser 50 pairs – Carton 400 pairs

FEELTOUCH 300PP SGA1SW

Хирургически ръкавици от естествен латекс, вътрешно покритие от царевично нишесте (USP), което се използва като лубрикант и за по-лесно поставяне.

Feeltouch предлага отлична комбинация от защита и комфорт, тъй като се вписва по ръката перфектно, благодарение на своята еластичност.

- Материал – естествен латекс
- AQL 1,0
- Леко напудрени
- Естествен цвят
- Микро текстурирана повърхност
- Анатомична форма (L / R) с извити пръсти
- Стерилизирани чрез гама-облъчване
- Медицинско изделие Class IIa 93/42/ЕИО
- Тест: EN 455 I-II-III, ASTM F 1671
- Размери 5,5-6-6,5-7-7,5-8-8,5-9
- Дозатор 50 чифта - кашон 400 чифта